


The PD Prescription

Ryan Blas BSN, RN

Clinical Cons. - Central Division



Disclaimer

- This presentation has been developed by Fresenius Medical Care North America.
- It is intended to provide pertinent data to assist healthcare professionals in forming their own conclusions and making decisions.
- It is not intended to replace the judgement or experience of the attending physician or other medical professionals.
- The treatment prescription is the sole responsibility of the attending physician.

Disclosure


I have nothing to disclose

PD Calculator Intended Use

The screenshot displays the '3. Physician Modeling' interface. It is divided into two main sections: 'Daytime' (marked with a sun icon) and 'Nighttime' (marked with a moon icon). Each section contains three input fields: 'Desired Fill Volume (L)', 'Desired Number Of Day/Night Exchanges', and 'Desired Time Per Exchange' (with a dropdown menu). Below these sections are two status indicators: 'total volume' (with a water drop icon) and 'total time' (with a clock icon). At the bottom, there is a gauge labeled 'Est. Total Weekly Kt/V' with a scale from 0 to 6 and a needle pointing to approximately 1.5.

- This course uses a PD Calculator as an example how modeling can facilitate finding individualized PD prescriptions.
- The online PD Calculator is intended for use by clinicians for modeling the dialysis dose (Kt/V) for different peritoneal dialysis prescriptions options. Modeling a PD dose is based on generalized formulas and assumptions derived from patient populations. The output of a modeled prescription is limited in its accuracy and cannot account for the variability seen in individual patients. It is essential that the physician adjusts the prescription according to the individual patient's clinical parameters to ensure the adequacy of the PD prescription.
- The calculator is not intended to replace the judgment or experience of the attending physician. The peritoneal dialysis treatment prescription is the sole responsibility of the attending physician.
- The PD Calculator is not intended to be used for pediatric patients or amputees.

Course Objectives



- ✓ Discuss adequacy of peritoneal dialysis
- ✓ Create an initial PD prescription using the PD Calculator
- ✓ Describe when and how to adjust and intensify the PD prescription
- ✓ Review sodium sieving and sodium removal in PD

Adequate Peritoneal Dialysis



Small solute clearance targets are met



Fluid balance, BP control, cardiac markers acceptable




Biochemical markers within acceptable ranges



Maximized quality of life

Factors That Affect Peritoneal Clearance



Transport Type

- Patient-specific
- Should be considered when determining appropriate dwell times for clearance and UF



Dialysate Volume

- More volume results in more clearance
- Can be adjusted by increasing # of exchanges or dwell volume



Dwell Time

- Should be based on transport type and clinical needs
- Avoid inappropriately short or long dwells
- Can be affected by catheter complications



Ultrafiltration Volume

- Difficult to predict accurately
- Individualizing dwell times for maximal UF is critical
- Adjust dextrose strength, consider icodextrin


Initial PD Prescriptions are Often Generic

... but your patients are not!!



Individualizing prescriptions is easy and a key to success in PD

Individualizing the PD Prescription



Case Study: New PD Patient

- ▶ Mr. R is an 82-year-old retired male new to dialysis
- ▶ Has some residual renal function with a GFR 6 mL/min*
- ▶ Height: 183 cm, Weight: 88 kg
- ▶ Transport type unknown
- ▶ Wants to start with manual PD

* For a patient of this size, a GFR ~6 mL/min equals a renal Kt/V ~ 1.34

Online PD Calculator

1. Patient Data

Age

Gender
Please Select ▼

Height
 cm in

Weight
 kg lb

Transport Status ⓘ
Please Select ▼

Residual Renal Function ⓘ
 K_{renal}/V mL/min


1. Enter basic patient data

- Age, gender, height, weight
- Transport status
- RRF

2. Estimated Prescription Data

Max. Fill Volume (L) ⓘ
--

Min. Number Of Exchanges (per day) ⓘ
--

Peak Time UF with 1.5% Glucose (hrs) ⓘ

Peak time shown in green. Time to negative UF shown in red.

Min. Total Daily Volume (L) ⓘ
--

BSA: --

Urea Distribution Volume: --

Renal Weekly Kt/V: --

2. Review estimated data

- Provides general guidance
- Based on modeling assumptions: max fill volume used, 4-hour dwells, continuous therapy, 1 L UF

3. Physician Modeling

Daytime ☀

Desired Fill Volume (L)

Desired Number Of Day Exchanges



Desired Time Per Exchange
Select ▼


Nighttime 🌙

Desired Fill Volume (L)

Desired Number Of Night Exchanges

Desired Time Per Exchange
Select ▼

 -- 
total volume total time

Est. Total Weekly Kt/V ⓘ


3. Predict outcomes

- Enter parameters and review modeling results for possible prescription options

4 Steps to Determine the PD Prescription

1 Enter patient-specific characteristics

2 Determine appropriate dwell volume

3 Consider clinical needs for approximating dwell times

4 Model # of exchanges and total volume required for adequate clearance

Determining the Initial Prescription

1 Patient Specifics

- Mr. R is an **82-year-old** retired **male** new to dialysis
- Height: **183 cm**, Weight: **88 kg**
- Transport type **unknown** – **assume average**
- Has decent residual renal function with a **GFR 6 mL/min**



1. Patient Data

Age

Gender

Height

 cm
 in

Weight

 kg
 lb

Transport Status ⓘ

Residual Renal Function ⓘ

 K_{renal}t/V
 mL/min

Determining the Initial Prescription

1 Patient Specifics

2 Dwell Volume

3

4

Dwell Volume

- Estimations for maximum:
 - Based on patient size
 - Includes UF
- Rule of Thumb:
 - 1.5 L/m² of BSA
 - Maximum of 3 L
- Initial volume for Mr. R: 1.5 – 2.0 L



2. Estimated Prescription Data


Max. Fill Volume (L) ⓘ

3 L * Initial volume is often chosen significantly lower and then increased over time

Min. Number Of Exchanges (per day) ⓘ

2

Peak Time UF with 1.5% Glucose (hrs) ⓘ



Peak time shown in green. Time to negative UF shown in red.

Min. Total Daily Volume (L) ⓘ

5.8

BSA: **2.1 m²**

Urea Distribution Volume: **45 L**

Renal Weekly Kt/V: **1.34**

Determining the Initial Prescription



1 Patient Specifics

2 Dwell Volume

3 Dwell Time

4

Dwell Time

- Peak time for UF:
 - Depends on transport type and PD solution used
- Appropriate time for Kt/V:
 - Correlates with peak time for UF
- Avoid:
 - Overly long dwells that result in reabsorption
 - Inefficient short cycles that result in inadequate UF, Kt/V, and sodium removal

2. Estimated Prescription Data


Max. Fill Volume (L) ⓘ

3 L

Min. Number Of Exchanges (per day) ⓘ

2

Peak Time UF with 1.5% Glucose (hrs) ⓘ



Peak time shown in green. Time to negative UF shown in red.

Min. Total Daily Volume (L) ⓘ

5.8

BSA: 2.1 m²

Urea Distribution Volume: 45 L

Renal Weekly Kt/V: 1.34

Determining the Initial Prescription

1 Patient Specifics

For many new PD patients with good RRF:

2 Dwell Volume

▶ Achieving adequate clearance is not difficult

BUT

3 Dwell Time


▶ Many initiate PD fluid overloaded

4


Consider a patient's UF requirements from the beginning

Fluid Status of Incident PD Patients

n = 1054



Dwell Time – Optimizing Ultrafiltration



Resource: Ultrafiltration Profiles


1.5% Dextrose

2.5% Dextrose


4.25% Dextrose

7.5% Icodextrin


Low Transporter




Low-Average Transporter



High-Average Transporter




High Transporter



Dwell Time – Optimizing Clearance



Clearance Profile for Average Transporter



Choosing a time around the peak time for UF seems a good approximation to ensure effective ultrafiltration and clearance





Resource: Clearance Profiles

1.5% Dextrose

2.5% Dextrose

4.25% Dextrose

7.5% Icodextrin



Determining the Initial Dwell Time



Dwell Time

1 Patient Specifics

- Mr R:
 - Good residual renal function

2 Dwell Volume

- Considerations:
 - UF peaks at 3 hours with 1.5% glc
 - Dwell times of 2 to 5 hr seem feasible
 - Only needs a few exchanges

3 Dwell Time

- Options:
 - Starting with a few 3 h manual exchanges with dry overnight
 - CAPD only possible with >1.5% glc

4

2. Estimated Prescription Data


Max. Fill Volume (L) ⓘ

3 L

Min. Number Of Exchanges (per day) ⓘ

2

Peak Time UF with 1.5% Glucose (hrs) ⓘ



Peak time shown in green. Time to negative UF shown in red.

Min. Total Daily Volume (L) ⓘ

5.8

BSA: 2.1 m²

Urea Distribution Volume: 45 L

Renal Weekly Kt/V: 1.34



Determining the Initial Prescription



2. Estimated Prescription Data

Max. Fill Volume (L) ⓘ
3 L

Min. Number Of Exchanges (per day) ⓘ
2

Peak Time UF with 1.5% Glucose (hrs) ⓘ

0 2 4 6 8 10 12

Peak time shown in green. Time to negative UF shown in red.

Min. Total Daily Volume (L) ⓘ
5.8

BSA: **2.1 m²**

Urea Distribution Volume: **45 L**

Renal Weekly Kt/V: **1.34**

3. Physician Modeling

Desired Fill Volume (L)

Desired Number Of Exchanges (per day)

Desired Time Per Exchange

Please Select ▼

-- --

total volume total time

Est. Total Weekly Kt/V ⓘ

--

Determining the Initial Prescription

1 Patient Specifics

2 Dwell Volume

3 Dwell Time



4 Model Options

3. Physician Modeling


Desired Fill Volume (L)

Desired Number Of Exchanges (per day)

Desired Time Per Exchange
 Please Select ▼

 -- 
 total volume total time

Est. Total Weekly Kt/V ⓘ



--

Modality Input ⓘ

Simple Day/Night



3. Physician Modeling

Daytime ☀


Desired Fill Volume (L)

Desired Number Of Day Exchanges

Desired Time Per Exchange
 Select ▼

 -- 
 total volume total time

Est. Total Weekly Kt/V ⓘ





--

Nighttime 🌙


Desired Fill Volume (L)

Desired Number Of Night Exchanges

Desired Time Per Exchange
 Select ▼

 -- 
 total volume total time


Est. Total Weekly Kt/V ⓘ



--

Modality Input ⓘ

Simple Day/Night





Example Initial Prescriptions

- 1 Patient Specifics
- 2 Dwell Volume
- 3 Dwell Time
- 4 Model Options

2. Estimated Prescription Data

Max. Fill Volume (L) **3 L**

Min. Number Of Exchanges (per day) **2**

Peak Time UF with 1.5% Glucose (hrs) **2**

0 2 4 6 8 10 12
Peak time shown in green. Time to negative UF shown in red.

Min. Total Daily Volume (L) **5.8**

BSA: **2.1 m²**

Urea Distribution Volume: **45 L**

Renal Weekly Kt/V: **1.34**

3. Physician Modeling

Daytime ☀	Nighttime 🌙
Desired Fill Volume (L) <input type="text"/>	Desired Fill Volume (L) 2
Desired Number Of Day Exchanges <input type="text"/>	Desired Number Of Night Exchanges 3
Desired Time Per Exchange Select	Desired Time Per Exchange 3 hours
6 L total volume	9 hours total time
Est. Total Weekly Kt/V 2.1	

3. Physician Modeling

Daytime ☀	Nighttime 🌙
Desired Fill Volume (L) 1.5	Desired Fill Volume (L) <input type="text"/>
Desired Number Of Day Exchanges 3	Desired Number Of Night Exchanges <input type="text"/>
Desired Time Per Exchange 5 hours	Desired Time Per Exchange Select
4.5 L total volume	15 hours total time
Est. Total Weekly Kt/V 2.1	

Example Initial Prescriptions

1

Patient
Specifics

Cycler Prescription

- Nightly Intermittent PD (NIPD)
 - 2 L fill volume
 - 3 exchanges per night
 - 3-hours per exchange
- Time on dialysis: 9-hr on cycler during the night
- Estimated Kt/V: 2.1
- Considerations:
 - Each exchange targeted for peak UF time with 1.5% dextrose
 - Adequate Kt/V
 - Free day

2

Dwell
Volume

3

Dwell
Time

4


Model
Options

Manual Prescription

- Manuals with dry overnight
 - 1.5 L fill volume
 - 3 exchanges per day
 - 5-hour exchanges
- Time on dialysis: 15 hours; i.e. from breakfast to dinner
- Estimated Kt/V: 2.1
- Considerations:
 - Low volume to provide comfort during day
 - Exchanges somewhat longer to adapt to time between breakfast and dinner
 - Dwell time still appropriate for 1.5% glucose
 - Adequate Kt/V
 - Free night, no need for cycler



Evaluating the Initial PD Prescription



**FRESENIUS
MEDICAL CARE**

4-8 Weeks After PD Start

Assess

Peritoneal transport status

Clearance & RRF

UF & Fluid Status

Clinical status

Intraperitoneal Pressure

Modify

Dwell times

Fill volumes

Glucose concentrations

PD Modality

Recommending Testing Timelines

Small Solute Clearance

- First month on dialysis
- At least every 4 months after
- Every 2 months in patients with residual renal function
- When clinically indicated

Peritoneal Membrane Status

- 4-8 weeks after PD initiation
- Changes in clearance or ultrafiltration
- Volume overload or increased need for hypertonic dialysate solutions
- Clinical symptoms of uremia or worsening hypertension


▶ K/DOQI . *Am J Kidney Dis.* 2006;48 Suppl 1:S91-S97

▶ Lo W-K, et al. *Perit Dial Int.* 2006;26(5):520-522

Determining Intraoperative Pressure

IPP > 18 cmH₂O becomes clinically symptomatic.
Normal IPP values range from 10-15 cmH₂O.

1. Have the patient lie in a perfectly horizontal position after instilling the dialysate
2. Set the zero level of the ruler to the mid-axillary line
3. Ensure infusion line is closed
4. Hang the drainage bag on a stand
5. Open the drain line and ensure the extension set is unclamped; the drain line will fill with effluent
6. Wait until the level of fluid in the drain line is stabilized with the patient breathing normally
7. Record IPP_{insp} and IPP_{exp} and average the two values
8. After measurements, lower the drain bag and drain the patient
9. Record drain volume



► Durand P-Y. *Perit Dial Int.* 2005;25(4):333-337

► Mathieu B, Pillot A. *Adv Perit Dial.* 1994;10:59-62

Case Study: Post-PET Adjustments



- ▶ Mr. R has been dialyzing for 6 weeks
- ▶ Adequacy Results:
 - **Renal Kt/V: 1.1**
 - **Peritoneal Kt/V: 0.8**
 - **Total Kt/V: 1.9**
- ▶ Height: 183 cm, Weight: 88 kg
- ▶ PET results: **High Average**
- ▶ Would consider to learning the cyclor

Case Study: Post-PET Adjustments

1 Patient Specifics

2

- Mr. R is an 82-year-old retired male
- Height: 183 cm, Weight: 88 kg
- Transport type: **High-Average**
- Renal **Kt/V: 1.1**

3

4



1. Patient Data

Age

82

Gender

Male

Height

183

 cm
 in

Weight

88

 kg
 lb

Transport Status ⓘ

High Average

Residual Renal Function ⓘ

1.1

 K_{renal}t/V
 mL/min


How does the time to peak UF change?

1 Patient Specifics

2 Dwell Volume

3 Dwell Time

4 Model Options



- A. Stays the same
- B. Gets shorter
- C. Gets longer




How does the time to peak UF change?

1 Patient Specifics

2 Dwell Volume

3 Dwell Time

4 Model Options



Should his current prescription be changed?



A. Yes
B. No

1 Patient Specifics

2 Dwell Volume

3 Dwell Time

4 Model Options

2. Estimated Prescription Data

Max. Fill Volume (L) **3 L**

Min. Number Of Exchanges (per day) **2**

Peak Time UF with 1.5% Glucose (hrs) **6**

Peak time shown in green. Time to negative UF shown in red.

Min. Total Daily Volume (L) **6.9**

BSA: **2.1 m²**

Urea Distribution Volume: **45 L**

Renal Weekly Kt/V: **1.1**

3. Physician Modeling


Daytime ☀	Nighttime 🌙
Desired Fill Volume (L) 1.5	Desired Fill Volume (L) <input type="text"/>
Desired Number Of Day Exchanges 3	Desired Number Of Night Exchanges <input type="text"/>
Desired Time Per Exchange 5 hours	Desired Time Per Exchange Select

4.5 L
total volume

15 hours
total time

Est. Total Weekly Kt/V **1.9**

General Modality Considerations Based on Transport Type, RRF and Body Size



Possible Adjusted Prescription



1 Patient Specifics

APD (NIPD)

- Increased fill volume
 - Since catheter healed in and supine position
- Shortened exchanges to 2.5 h
 - Accounts for transport status
 - Adjusts for max UF per 1.5% glucose exchange
 - Dry day
- Time on dialysis: 7.5 hr
- Estimated Kt/V: 2.2


2 Dwell Volume


3 Dwell Time

4 Model Options


3. Physician Modeling

Daytime ☀	Nighttime 🌙
Desired Fill Volume (L) <input type="text"/>	Desired Fill Volume (L) <input type="text" value="2.5"/>
Desired Number Of Day Exchanges <input type="text"/>	Desired Number Of Night Exchanges <input type="text" value="3"/>
Desired Time Per Exchange <input type="text" value="Select"/>	Desired Time Per Exchange <input type="text" value="2.5 hours"/>

 **7.5 L**
total volume


 **7.5 hours**
total time

Est. Total Weekly Kt/V ⓘ



2.2

Intensifying the PD Prescription



When to Intensify Dialysis

Volume overload

Loss of RRF

Hyperkalemia

Hyperphosphatemia

Metabolic acidosis

Anemia

Uremic neuropathy

Uremic pericarditis

Pruritus

Nausea or vomiting

Sleep disturbances
(restless leg syndrome)

Intensifying the PD Prescription


Clearance

- ↑ # of exchanges
- Adjust dwell time
- ↑ Fill volume
- Change modality

Ultrafiltration

- Adjust dwell time
- ↑ Osmotic strength
- ↓ IPP
- Dietary restrictions

Case Study: Loss of RRF



- ▶ Mr. R has been dialyzing for 15 months
- ▶ Adequacy Results:
 - **Renal Kt/V: 0**
 - **Peritoneal Kt/V: 1.1**
 - **Total Kt/V: 1.1***
- ▶ Height: 183 cm, Weight: 88 kg
- ▶ PET results: **High Average**

* Obviously, this prescription should have been adjusted earlier. This is a didactic example to demonstrate how much Kt/V may decrease if prescription is not regularly adjusted.

Case Study: Lost RRF

1 Patient Specifics

2

3

4

- Mr. R is an 83-year-old, retired
- Height: 183 cm, Weight: 88 kg
- Transport type: **High-Average**
- Renal **Kt/V: 0**



1. Patient Data

Age

82

Gender

Male

Height

183

 cm
 in

Weight

88

 kg
 lb

Transport Status ⓘ

High Average

Residual Renal Function ⓘ

0

 K_{renal}t/V
 mL/min

What prescription would you pick?



1 Patient Specifics

A. CCPD with 4 x 2.8 L, 2.5-hr dwells with 1.5% glc; last fill of 2.3 L with 4.25% or icodextrin

2 Dwell Volume

B. PD Plus™ with 4 x 2.5 L, 2-hr dwells with 1.5% glc; 2-day exchanges of 2L with 2.5% glc

3 Dwell Time

C. CCPD with 6 x 2.5 L, 1.5-hr dwells with 1.5% glc; last fill of 2.0 L with 4.25% glc or icodextrin

4 Model Options

D. NIPD with 8 x 3 L, 1-hr dwells with 1.5% glc; dry day

Possible PD Plus Prescription

1 Patient Specifics

2 Dwell Volume

3 Dwell Time

4 Model Options



- Mr. R wants to try PD Plus
- Since retired and home he does not mind the extra day exchange
- He likes the shorter night and low fill volumes

2. Estimated Prescription Data

Max. Fill Volume (L) ⓘ
3 L

Min. Number Of Exchanges (per day) ⓘ
6

Peak Time UF with 1.5% Glucose (hrs) ⓘ

0 2 4 6 8 10 12
Peak time shown in green. Time to negative UF shown in red.

Min. Total Daily Volume (L) ⓘ
15.4

BSA: **2.1 m²**

Urea Distribution Volume: **45 L**

Renal Weekly Kt/V: **0**

3. Physician Modeling

Daytime ☀	Nighttime 🌙
Desired Fill Volume (L) <input type="text" value="2"/>	Desired Fill Volume (L) <input type="text" value="2.5"/>
Desired Number Of Day Exchanges <input type="text" value="2"/>	Desired Number Of Night Exchanges <input type="text" value="4"/>
Desired Time Per Exchange <input type="text" value="8 hours"/>	Desired Time Per Exchange <input type="text" value="2 hours"/>


14 L
total volume

24 hours
total time

Est. Total Weekly Kt/V ⓘ

1.9

Case Study: Hospitalization for Fluid Overload




- ▶ Mr. R has been feeling increasingly ill and presents with a pronounced increase in weight
- ▶ Admitted to the hospital for complications related to severe fluid overload
- ▶ Height: 183 cm, Weight: **93 kg (dry weight 88)**
 - Resident nephrologist prescribes 6 x 2 L, 1-hr dwells with 4.25% dextrose to pull fluid off
 - After 4 exchanges, Mr. R complains of severe thirst, nausea, and muscle spasms

What happened?



- A. Fluid was removed too fast and patient became hypovolemic
- B. The dwell times were too short and patient became hypernatremic
- C. The dwell times were too long and he remained hypervolemic
- D. None of the above.



Sodium Sieving in Peritoneal Dialysis



- Due to free water transport at the beginning of a dwell:
 - Dialysate sodium is diluted
 - Plasma sodium is “sieved” and becomes more concentrated
 - Relevant especially early in the dwell
- Note: sodium sieving looks at the *concentration* of sodium not at the *amount* removed

Sodium Removal in PD

- Mainly depends on 2 variables:
 - Sodium sieving and Ultrafiltration
- Sieving peaks in the first 1-2 hours and then steadily improves
- This seems to indicate the longer the dwell the better the sodium removal
- However sodium removal depends most strongly on UF
 - Total sodium removal declines in long dwells after peak UF is reached
- Sieving impairs sodium removal in too short cycles, reabsorption in too long cycles




Sodium Removal



- ▶ Sodium removal correlates heavily with ultrafiltration volume
 - No significant concentration gradient for sodium
 - Sodium removal thus largely dependent on convective transport
 - Convective transport depends on UF rate
 - Sodium removal peaks around the same time peak ultrafiltration is reached with each dextrose concentration
- ▶ **Bottom line:** Maximizing prescription for time to peak UF also maximizes sodium removal

▶ Wang T, et al. *Kidney Int.* 1997;52(6):1609-1616

Correlation of Sodium Removal with UF




Case Study: Hospitalization for Fluid Overload




- The dwell time of 1-hr with a 4.25% glucose solution resulted in marked sodium sieving as indicated by the red crosses
- The repeated short cycles with high dextrose solution lead to hypernatremia
- A change in the regimen to 4-hr dwells achieved adequate sodium and fluid removal and resolved the hypernatremia
- Mr. R was discharged after 2 days

Dietary Sodium Restriction

Reinforce importance of dietary sodium intake and fluid restrictions



Clinical Takeaways



- ✓ Don't assume 4 x 2L exchanges is ideal for all patients
- ✓ Customize prescription to clinical needs and patient transport type
- ✓ Assess more than Kt/V when determining adequacy
- ▶ Optimize dwell time for peak UF which is also best for sodium removal
- ✓ Avoid very short dwell times that may result in sodium sieving
- ▶ Avoid dwell times that are too long and result in fluid reabsorption and negative UF

Questions?

